



قائمة الاسئلة

تقويم الاسنان(2)-الثالث -فني الاسنان- كلية طب وجراحة الفم والاسنان- درجة الامتحان (30)

د. صالح البقال

- 1) Reverse loop canine retractor is indicated when :
 - 1) - The canine is placed palatally
 - 2) + The canine is placed in the line of the arch and has to be just distalized
 - 3) - The functional depth of the sulcus is Less as in the lower arch
 - 4) - The canine is placed labially or high in the sulcus
- 2) It is indicated when the minimal canine retraction of 1 to 2 mm is required :
 - 1) - Reverse Loop Canine Retractor
 - 2) - Buccal Canine Retractor
 - 3) - Palatal Canine Retractor
 - 4) + "U" Loop Canine Retractor
- 3) We use 0.6 mm round S.S wire for it:
 - 1) - Helical Canine Retractor
 - 2) - "U" Loop Canine Retractor
 - 3) - Buccal Canine Retractor
 - 4) + Palatal Canine Retractor
- 4) To Fabricate a Supported Buccal Canine Retractor use :
 - 1) - Wire of 0.7 mm
 - 2) - Wire of 0.9 mm
 - 3) + Wire of 0.5 mm + Tube
 - 4) - Wire of 0.7 mm + Tube
- 5) To Activate the Palatal Canine Retractor, We :
 - 1) - Close the Coil 2-3 mm
 - 2) + Open the Coil 2-3 mm
 - 3) - Don't need Activation
 - 4) - Bending the Guide arm
- 6) When we Activate the jackscrew, per quarter-turn the screw produces :
 - 1) - 2 - 2.5 mm movement
 - 2) - 3 - 4 mm movement
 - 3) + 0.2 - 0.25 mm movement
 - 4) - 0.7 - 0.75 mm movement
- 7) The manufactured screw is designed to be opened by one-quarter turn, that is less than the total width of periodontal ligament that are :
 - 1) + about 0.15 mm to 0.35 mm
 - 2) - about 1.7 mm to 1.75 mm
 - 3) - about 0.35 mm to 0.50 mm
 - 4) - about 2.15 mm to 2.75 mm
- 8) The three Dimensional expansion screw is usually used for :
 - 1) - Closure of extraction spaces
 - 2) + Class III cases
 - 3) - Split activator appliance
 - 4) - Lateral and sagittal anterior expansion of the mandible arch
- 9) The Telescopic Screw is usefully used for :
 - 1) - Class II div I cases
 - 2) - Bilateral symmetrical expansion of the maxillary arch
 - 3) - Reduction over jet



- 4) ☒ Moving a single tooth in the labial direction
- 10) The Sagittal expansion Appliance is :
- 1) ☐ Bi-Lateral Symmetrical Expansion
 - 2) ☐ Radial Expansion
 - 3) ☒ Anterior/posterior expansion
 - 4) ☐ Not all of these
- 11) Anterior Crossbite is a malocclusion occurring in which plane :
- 1) ☐ Vertical Plane
 - 2) ☐ Transverse Plane
 - 3) ☒ Sagittal Plane
 - 4) ☐ Horizontal Plane
- 12) All of the following are Ways for Expansion, EXCEPT :
- 1) ☐ Ni-Ti palatal expander
 - 2) ☐ Screws
 - 3) ☐ Quad Helix Appliance
 - 4) ☒ Nance Appliance
- 13) Quad Helix Appliance is a:
- 1) ☒ Expansion Appliance
 - 2) ☐ Space Maintainer Appliance
 - 3) ☐ Habit Breaker Appliance
 - 4) ☐ Passive Appliance
- 14) Traction Screws are Used for :
- 1) ☒ Closure of Extraction spaces
 - 2) ☐ Expansion Dental Arch
 - 3) ☐ Correction a Single Cross bite
 - 4) ☐ Activator Appliances
- 15) Bertoni & Steiner Screws are a Type of :
- 1) ☐ Micro Screws
 - 2) ☐ S.S Wires
 - 3) ☒ Three-Dimensional Expansion Screws
 - 4) ☐ Fan Expansion Screws
- 16) Hyrax Appliance is considerate a :
- 1) ☒ Tooth borne appliance
 - 2) ☐ Tooth tissue borne appliances
 - 3) ☐ Removable Appliance
 - 4) ☐ Habit Breaker Appliance
- 17) The Appliance that can only be Used For Slow Expansion is :
- 1) ☒ Coffin Spring
 - 2) ☐ Hass Expander
 - 3) ☐ Hyrax Screw
 - 4) ☐ Fan RME type
- 18) The main Indication of Anterior Bite Plane is :
- 1) ☐ Class II cases
 - 2) ☐ Class III cases
 - 3) ☐ Open bite cases
 - 4) ☒ Deep bite cases
- 19) The Indications of Posterior Bite Plane is all of the following Except :
- 1) ☐ Cross Bite
 - 2) ☐ TMJ pain dysfunction cases



- 3) ☒ Deep bite cases
4) ☐ Slight Open bite cases
- 20) Reduction of overbite is done by :
1) ☒ Anterior bite plane
2) ☐ Posterior bite plane
3) ☐ Both anterior & posterior bite plane
4) ☐ Labial bow with Adam's clasp & Z-spring
- 21) For cases of Thumb Sucking, we can use all of these, EXCEPT:
1) ☐ Boxing gloves
2) ☒ Lingual Bar
3) ☐ Thumb crib
4) ☐ Finger cab
- 22) What is Bluegrass Appliance ?
1) ☐ Space Maintainer
2) ☐ Expansion Appliance
3) ☒ Habit Breaker
4) ☐ Active Appliance
- 23) Oral Screen is used to correct the :
1) ☒ Mouth Breathing
2) ☐ Bruxism
3) ☐ Lip Biting Habit
4) ☐ Tongue thrusting
- 24) The Space Maintainer is considered a type of :
1) ☐ Retainers Orthodontics
2) ☐ Myofunctional Orthodontics
3) ☐ Habit Breakers
4) ☒ Preventive Orthodontics
- 25) The Best Space Maintainers a:
1) ☐ Band & loop
2) ☒ Primary Tooth
3) ☐ Distal Shoe
4) ☐ Nance
- 26) When the deciduous two centrals (A,A) is lost or extracted at the age of 7 years, what is the suggested space Maintainers?
1) ☐ Pado Partial
2) ☐ Nance
3) ☐ Band & loop
4) ☒ No need
- 27) The "Leeway Space" in the Maxilla is usually about:
1) ☒ 1 to 2 mm
2) ☐ 2 to 3 mm
3) ☐ 3 to 4 mm
4) ☐ 4 to 5 mm
- 28) One of the functions of the space maintainers is to prevent tongue thrust.
1) ☒ TRUE.
2) ☐ FALSE.
- 29) What is TPA?
1) ☐ Treatment planning approach
2) ☐ Transverse plane articulation



- 3) ☒ Transpalatal arch
4) ☐ Triangular palatal activation
- 30) The correct dentition in 5 years old child in the upper jaw is :
- 1) ☒ ABCDE
2) ☐ 12CDE6
3) ☐ 123DE6
4) ☐ 12C456