



قائمة الاسئلة 2025-05-03 06:51

زراعات تعويضية-الثالث-تقني اسنان- كلية طب وجراحة الفم والاسنان - درجة الامتحان(25)

Dr.Maha Al-khorasani

- 1) are instruments that provide predictable retention of abutments for extr-aoral preparation and polishing procedures
  - 1) - Replicas
  - 2) ☒ Abutment Holders
  - 3) - Cylinders
  - 4) - none
- 2) CAD/CAM Abutments may
  - 1) - Milled Zirconia
  - 2) ☒ all
  - 3) - Titanium
  - 4) - none
- 3) UCLA Abutments
  - 1) - are ideal for tissue supported overdentures on two to four implants
  - 2) - all
  - 3) - nothing
  - 4) ☒ developed to allow a direct connection between implants and implant restorations
- 4) Bar attachment is
  - 1) ☒ An Abutment used with removable prosthesis
  - 2) - An Abutment used with fixed prosthesis
  - 3) - A Healing Abutments
  - 4) - nothing
- 5) Healing abutment
  - 1) - is designed to be part of Emergence Profile to guide soft tissue healing after implant placement
  - 2) - is available for all implant restorative platforms
  - 3) ☒ all
  - 4) - is designed with collar heights between 2 and 8 mm and 3.4, 4.1, 5, 6, and 7.5 mm diameters
- 6) is the portion of the implant that supports or retains a prosthesis or implant superstructure
  - 1) - Fixture of Implant
  - 2) - Crest Module
  - 3) ☒ Abutment of implant
  - 4) - none
- 7) Wide Diameter implants are usually
  - 1) - 3.5 - 4.5 mm diameter
  - 2) ☒ > 4.5 mm in diameter
  - 3) - < 3.5 mm
  - 4) - none
- 8) Combination implants
  - 1) ☒ usually are pressed into position and have a macro body design similar to a thread form
  - 2) - are non-threaded
  - 3) - All
  - 4) - none
- 9) The implant body is consist of
  - 1) - crest module
  - 2) - body
  - 3) - apex



- 4) ☒ All
- 10) Transosteal implant subdivided into
- 1) ☒ Staple, Single pin and Multiple pin
  - 2) ☐ Root form, Blade (plate) form and Ramus frame
  - 3) ☐ all
  - 4) ☐ none
- 11) Contraindications for implant treatment
- 1) ☒ all
  - 2) ☐ Untreated intraoral pathology or malignancy
  - 3) ☐ Untreated periodontal disease
  - 4) ☐ Uncontrolled drug or alcohol use
- 12) Disadvantages of Implant-Supported Protheses
- 1) ☐ Maintain facial esthetics
  - 2) ☒ none
  - 3) ☐ Restore and maintain occlusal vertical dimension
  - 4) ☐ Improve phonetics
- 13) exact replicas of implant abutments and are manufactured to the same high degree of accuracy as implant analogs
- 1) ☐ Implant Analogs
  - 2) ☒ Abutment Analogs
  - 3) ☐ Abutment Holders
  - 4) ☐ none
- 14) LOCATOR Abutments are ideal for
- 1) ☐ all
  - 2) ☐ crown and bridges
  - 3) ☒ tissue supported overdentures
  - 4) ☐ none
- 15) implant abutment can be classified according to their shape into
- 1) ☐ All
  - 2) ☐ none
  - 3) ☐ Ready made or Castable
  - 4) ☒ Straight or Angulated
- 16) Types of Healing Abutment
- 1) ☐ Stock
  - 2) ☒ All
  - 3) ☐ CAD/CAM
  - 4) ☐ Custom
- 17) Anti-rotation feature of platform are
- 1) ☐ External hex
  - 2) ☒ All
  - 3) ☐ Internal hex
  - 4) ☐ none
- 18) it represents the transition zone from the implant body design to the transosteal region of the implant at the crest of the ridge
- 1) ☐ none
  - 2) ☐ Fixture of Implant
  - 3) ☐ platform
  - 4) ☒ Crest Module
- 19) Implants are generally available in lengths from



- 1) - none
  - 2) - about 10 mm to as much as 20 mm
  - 3) - about 2 mm to as much as 20 mm
  - 4) ☒ about 6 mm to as much as 20 mm
- 20) implant abutment may be classified according to the way of abutment connection
- 1) ☒ Screwed and Cemented
  - 2) - Straight and Angulated
  - 3) - Ready made and Castable
  - 4) - none
- 21) These abutments were manufactured and designed for use with screw-retained restorations
- 1) ☒ Low-Profile Conical Abutments
  - 2) - LOCATOR
  - 3) - All
  - 4) - none
- 22) CAD/CAM Abutments may be
- 1) - none
  - 2) - Titanium
  - 3) - Milled Zirconia
  - 4) ☒ All
- 23) is used for screw insertion and removal
- 1) ☒ Drivers and Placement Instruments
  - 2) - protector
  - 3) - analog
  - 4) - none
- 24) Options of treatment for tooth loss
- 1) - No treatment
  - 2) - A bridge
  - 3) - removable dentur
  - 4) ☒ All
- 25) it is a negative Effect of Removable Prostheses
- 1) ☒ The life span may be decreased
  - 2) - The life span may be increased
  - 3) - All
  - 4) - none